



Recombinant Human Zinc finger protein 654 (ZNF654)

Product Code	CSB-EP026922HU
Abbreviation	ZNF654
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZM8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLCNKEFLGG HIVRHAQAHQ KKGSFACVIC GRKFRNRGLM QKHLKNHVKK IQRQQIAAAQ QDDQEVTALE EINCSSSSIS FENGNSDSKD LEVETLTASS EGNKEVIPEH VAEFIEIPIS VPEDVIENVI ENGSPNNSLN NVFKPLTECG DDYEEEEDEE GDYEEDDYDL NQETSVIHKI NGTVCHPKDI YATDQEGNFK CPALGCVRIF KRIGFLNKHA MTVHPTDLNV RQTMKWSKG KCKFCQRQFE DSQHFIDHLN RHSYPNVYFC LHFNCNESFK LPFQLAQHTK SHRIFQAQCS FPECHELFED LPLLYEHEAQ HYLKTPESS AQPSETILWD VQTDSNPNQE KDSSSNEKQT ISLPVSTSKS RKESTEPKTC IESMEKKTDS LVQNGNERSD DTVSNISLID QKMPDIEPNS ENNCSSSDIV NGHSEIEQTP LVSSDPALKI DTNRIRTENG SILPSVVPQE HNTLPVSQAP SKPNLTSEHT SYGLILTKPY VRPLPPSYLD ERYLSMPKRR KFLTDRVDAC SDQDNVYKKS VKRLRCGKCL TTYCNAEALE AHLAQKQCQT LFGFDSDES A
Source	E.coli
Target Names	ZNF654
Protein Names	Recommended name: Zinc finger protein 654 Alternative name(s): Melanoma-associated antigen
Expression Region	1-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.