



Recombinant Human Cytokine receptor-like factor 3 (CRLF3)

Product Code	CSB-MP005984HU
Abbreviation	CRLF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUI8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRGAMELEPE LLLQEARENV EAAQSYRREL GHRLEGLREA RRQIKESASQ TRDVLKQHFN DLKGTGKLL DERLVTLLQE VDTIEQETIK PLDDCQKLIE HGVNTAEDLV REGEIAMLGG VGEENEKLWS FTKKASHIQL DSLPEVPLLV DVPCLSAQLD DSILNIVKDH IFKHGTVASR PPVQIEELIE KPGGIIVRWC KVDDDFTAQD YRLQFRKCTS NHFEDVYVGS ETEFIVLHID PNVDYQFRVC ARGDGRQEWS PWSVPQIGHS TLPHEWTAG FEGYSLSSRR NIALRNDSES SGVLYSRAPT YFCGQTLTFR VETVGQPD RR DSIGVCAEKQ DGYDSLQRDQ AVCISTNGAV FVNGKEMTNQ LPAVTSGSTV TFDIEAVTLG TTSNNEGGHF KLRVTISSNN REVVFDWLLD QSCGSLYFGC SFFYPGWKVL VF
Source	Mammalian cell
Target Names	CRLF3
Protein Names	Recommended name: Cytokine receptor-like factor 3 Alternative name(s): Cytokine receptor-like molecule 9 Short name= CREME-9 Cytokine receptor-related protein 4 Type I cytokine receptor-like factor p48
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.