



Recombinant Human Putative homeodomain transcription factor 2 (PHTF2)

Product Code	CSB-MP017945HU
Abbreviation	PHTF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3S3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASKVTDIV WYQKKEFLSV ATTAPGPQQV LPGYCQCSLK DQGLFIQCLI GAYDQQIWEK SVEQREIKFI KLGLRNKPKK TAHVKPDLID VDLVRGSAFA KAKPESPWTS LTRKGIVRVV FFPFFFRWWL QVTSKVIFFW LLVLYLLQVA AIVLFCSTSS PHSIPLTEVI GPIWLMLLLG TVHCQIVSTR TPKPPLSTGG KRRRKLKRAA HLEVHREGDG SSTTDNTQEG AVQNHGTSTS HSVGTVFRDL WHAFFLSGS KKAKNSIDKS TETDNGYVSL DGKKTVKSGE DGIQNHEPQC ETIRPEETA W NTGTLRNGPS KDTQRTITNV SDEVSSSEEGP ETGYSLRRHV DRTSEGVLRN RKSHHYKKHY PNEDAPKSGT SCSSRCSSSR QDSESARPES ETEDVLWEDL LHCAECHSSC TSETDVENHQ INPCVKKEYR DDPFHQSHLP WLHSSHPGLE KISAIWEGN DCKKADMSVL EISGMIMNRV NSHIPGIGYQ IFGNAVSLIL GLTPFVFRLS QATDLEQLTA HSAELYVIA FGSNEDVIVL SMVIISFVVR VSLVWIFFFL LCVAERTYKQ RLLFAKLFGH LTSARRARKS EVPHFRLKKV QNIKMWLSLR SYLKRRGPQR SVDVIVSSAF LLTISVVFIC CAQLLHVHEI FLDCHYNWEL VIWCISLTLF LLRFVTLGSE TSKKYSNTSI LLTEQINLYL KMEKKPNKKE ELTLVNNVLK LATKLLKELD SPFRLYGLTM NPLLYNITQV VILSAVSGVI SDLLGFNLKL WKIKS
Source	Mammalian cell
Target Names	PHTF2
Protein Names	Recommended name: Putative homeodomain transcription factor 2
Expression Region	1-785
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.