



Recombinant Human Zinc finger protein 433 (ZNF433)

Product Code	CSB-EP026747HU-B
Abbreviation	ZNF433
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N7K0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMFQDSVAFE DVAVTFTQEE WALLDPSQKN LCRDVMQETF RNLASIGKKW KPQNIYVEYE NLRRNLRIVG ERLFESKEGH QHGEILTQVP DDMLKKTITG VKSCCESSVYG EVGSAHSSLN RHIRDDTGHK AYEYQEYQK PYKCKYCKKP FNCLSSVQTH ERAHSGRKLY VCEECGKTFI SHSNLQRHRI MHRGDGPYKC KFCGKALMFL SLYLIHKRTH TGEKPYQCKQ CGKAFSHSSS LRIHERHTG EKPYKNECG KAFHSSTCLH AHKRTHTGEK PYECKQCGKA FSSSHSFQIH ERTHTGEKPY ECKECGKAFK CPSSVRRHER THSRKKPYEC KHCGKVL Syl TSFQNHLMGMH TGEISHKCKI CGKAFYSPSS LQTHEKTHTG EKP YKCNQCG KAFNSSSSFR YHERHTHTGEK PYECKQCGKA FRSASLLQTH GRHTHTGEKPY ACKECGKPFNSFSFFQIHER MHREEKPYEC KGYGKTFSLP SLFHRHERTH TGGKTYECKQ CGRSFNCSSS FRYHGRHTHTG EKP YECKQCG KAFRSASQLQ IHGRHTHTGEK PYECKQCGKA FGSASHLQMH GRHTHTGEKPY ECKQCGKSFG CASRLQMHGR THTGEKPYKC KQCGKAFGCP SNLRRHGRTH TGEKPYKCNQ CGKVFRCSQ LQVHGRAHCI DTP
Source	E.coli
Target Names	ZNF433
Protein Names	Recommended name: Zinc finger protein 433
Expression Region	1-673
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.