



Recombinant Human Zinc finger protein 548 (ZNF548)

Product Code	CSB-EP026829HU-B
Abbreviation	ZNF548
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEK5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNLTEGRVVF EDVAIFYFSQE EWGHLDEAQR LLYRDVMLEN LALLSSLGSW HGAEDDEEAPS QQGFSVGVSE VTASKPCLSS QKVHPSETCG PPLKDILCLV EHNGIHPEQH IYICEAELFQ HPKQQIGENL SRGDDWIPSF GKNHRVHMAE EIFTCMEGWK DLPATSCLLQ HQGPQSEWKP YRDTEDREAF QTGQNDYKCS ECGKTFTCSY SFVEHQKIHT GERSYECNKC GKFFKYSANF MKHQTVHTSE RTYECRECGK SFMYNYRLMR HKRVHTGERP YECNTCGKFF RYSSTFVRHQ RVHTGERPYE CRECGKFFMD SSTLIKHQRV HTGERPYKCN DCGKFFRYIS TLIRHQRIHT GERPYECSVC GELFRYNSSL VKHWRNHTGE RPYKCSECGK SFRYHCRLIR HQRVHTGERP YECSECGKFF RYNSNLIKHW RNHTGERPYE CRECGKAFSH KHILVEHQKI HSGERPYECS ECQKAFIRKS HLVHHQKIHS EERLVCSMNV GNSLAKTPTS LNIRDFTMEK VYH
Source	E.coli
Target Names	ZNF548
Protein Names	Recommended name: Zinc finger protein 548
Expression Region	1-533
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.