



Recombinant Human Folliculin (FLCN)

Product Code	CSB-EP008711HU-B
Abbreviation	FLCN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NFG4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNAIVALCHFCELGHPRTLFCTEVLHAPLPQGDGNE SPGQGEQAE EEEEGGI QMNSRMRA HSPAEGASVESSSPGPKKSDMCEGCRSLAAGHPGYISHDKETS SIKYVSHQHP SHPQLFSI VRQACVRSLSCEVCPGREGPIFFGDEQHG FVSHTFFIKDSLARGFQRWYSIIT IMMDRI YLINSWPFLGKVRGIIDELQGKALKVFEAEQFGCPQRAQRMNTAFTPFLHQR NGNAARS LTSLTSDDNLWACLHTSF AWLLKACGSRLTEKLLEGAPTE DTLVQMEKLAGEA GVLLPGP WPGWPWGGTSCLLSWQESLREGNAALNQPRTSLREAHPPISV
Source	E.coli
Target Names	FLCN
Protein Names	Recommended name: Folliculin Alternative name(s): BHD skin lesion fibrofolliculoma protein Birt-Hogg-Dube syndrome protein
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is located within the Smith-Magenis syndrome region on chromosome 17. Mutations in this gene are associated with Birt-Hogg-Dube syndrome, which is characterized by fibrofolliculomas, renal tumors, lung cysts, and pneumothorax. Alternative splicing of this gene results in two transcript variants encoding different isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.