



Recombinant Human Zinc finger protein 502 (ZNF502)

Product Code	CSB-EP026795HU-B
Abbreviation	ZNF502
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TBZ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLNMQGAEER DIRRETCPGW VNKNKPALEQ DVCKIDSSGI VVKRFQEDEY QDSTFEEKYA CEGMKENSPR EIAESCLFQE GGFGRITFIH KEAPPEIISQ GYNFEKSLLL TSSLVTRLRV STEESLHQWE TSNIQTNDIS DQSKCPTLCT QKKSWKCNEC GKTFTQSSSL TQHQRTHTGE RPYTCEECEGK AFSRSSFLVQ HQRIHTGVKP YGCEQCGKTF RCRSFLTQHQ RIHTGEKPYK CNECGNSFRN HSHLTEHQRI HTGEKPYKCN RCGKAFNQNT HLIHHQRIHT GEKPYICSEC GSSFRKHSNL TQHQRHTGE KPHKCDECGK TFQTKANLSQ HQRIHSGEKP YKCKEKGKAF CQSPSLIKHQ RIHTGEKPYK CKECGKAFTQ STPLTKHQRI HTGERPYKCS ECGKAFIQSI CLIRHQRSHT GEKPYKCNEC GKGFNQNTCL TQHMRIHTGE KPYKCKEKGK AFAHSSSLTE HHRHTHTGEKL YKCSECEKTF RKYAHLSEHY RIHTGEKPYE CIECGKFFRH SSVLFRHQKL HSGD
Source	E.coli
Target Names	ZNF502
Protein Names	Recommended name: Zinc finger protein 502
Expression Region	1-544
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.