



# Recombinant Mouse Suppressor of cytokine signaling 7 (Socs7)

<b>Product Code</b>	CSB-EP022396MO-B
<b>Abbreviation</b>	Socs7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8VHQ2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVFRNVGRPP EEEDAEAARE PGPSELLCPR HRCALDPKAL PPGLALERTW GPVAGLEAQL AALGLGQPAG PGIKTAGGGC CPCPCPPQPP PPQPPPPAAA PQAGEDPTET SDALLVLEGL ESEAESLETN SCSEELSSP GRGGGGVGGR LLLQPPGPEL PPVPFPLQDL VPPGRLSRGE QQQQQPPPPP PPPGPLRPLA GPSRKGSFKI RLSRLFRTKS CNGGSGGGDG TGKRPSGDLA ASAASLTDMG GSAVRELDTG RKPRLTRTQS AFSPVSFSPL FTGETVSLVD VDISQRGLTS PHPTPPPPP RRSLSLLDDI SGTLPSTVLV APMGSSLQSF PLPPPPPPHA PDAFPRIAPI RASESLHSQP PQHLQCPLYR PDSSSFAASL RELEKCGWYW GPMNWEDAEM KLKKGKPDGSF LVRDSSDPY ILSLSFRSQG ITHHTRMEHY RGTFSWLWCHP KFEDRCQSVV EFIKRAIMHS KNGKFLYFLR SRVPGLPPTP VQLLYPVSRF SNVKSQHLR RFRIRQLVRI DHIPDLPLPK PLISYIRKFY YYDPQEEVYL SLKEAQLISK QKQVEVPST
<b>Source</b>	E.coli
<b>Target Names</b>	Socs7
<b>Protein Names</b>	Recommended name: Suppressor of cytokine signaling 7 Short name= SOCS-7
<b>Expression Region</b>	1-579
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.