



Recombinant Human Heterogeneous nuclear ribonucleoprotein L-like (HNRPLL)

Product Code	CSB-EP010630HU
Abbreviation	HNRPLL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WVV9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSSSSSPRE TYEEDREYES QAKRLKTEEG EIDYSAEEGE NRREATPRGG GDGGGGGRSF SQPEAGGSHH KVSVSPVVHV RGLCESVVEA DLVEALEKFG TICYVMMMPF KRQALVEFEN IDSAKECVTF AADEPVYIAG QQAFFNYSTS KRITRPGNTD DPSGGNKVLL LSIQNPLYPI TVDVLYTVCN PVGKVQRIVI FKRNGIQAMV EFESVLCAQK AKAALNGADI YAGCCTLKIE YARPTRLNVI RNDNDSWDYT KPYLGRRDRG KGRQRQAILG EHPSSFRHDG YGSHGPLLPL PSRYRMGSRD TPELVAYPLP QASSSYMHGG NPSGSVVMVS GLHQLKMNCS RVFNLFCLYG NIEKVKFMKT IPGTALVEMG DEYAVERA VAVT HLNNVKLFGK RLNVCVSKQH SVVPSQIFEL EDGTSSYKDF AMSKNNRFTS AGQASKNIIQ PPSCVLHYYN VPLCVTEETF TKLCNDHEVL TFIKYKVFDA KPSAKTLSGL LEWECKTDAV EALTALNHYQ IRVPNGSNPY TLKLCFSTSS HL
Source	E.coli
Target Names	HNRNPLL
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein L-like Short name= hnRNPLL Alternative name(s): Stromal RNA-regulating factor
Expression Region	1-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.