



Recombinant Bovine Dentin matrix acidic phosphoprotein 1 (DMP1)

| | |
|--------------------------|--|
| Product Code | CSB-YP006967BO |
| Abbreviation | DMP1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q95120 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | LPVA RYQNTESKSS EEWKGHLAQT PTPPLESSES SEESKLSSEE QANEDPSDST ESEEVGLDD QQHVHRPAGG LSRRGGSEGD NKDDDEDESG DDTFGDDGG PGPEERRSGG DSRLGSDSDS ADTTRSREDS TPQGDEGARD TTSESRLDR EDEGNSRPEG GDSTPDSSE EHWVGGGSEG DSSHGDGSEF DDEGMQSDDP GAYRSERGNS RISDAGLKST QSKGDDEEQA STQDSHESPA AAYPRRKFFR KSRLPEEDGR GELDDSRITIE VMSDSTENPD SKEAGLGQSR EHSKSESRQE SEENRSPEDS QDVQDPSSSES SQEVDLPSQE NSSSESQEEAL HESRGDNPDN ATSHSREHQA DSESSEEDVL DKPSDSESTS TEEQADSESH ESLRSSEESP ESTEEQNSSS QEGAQTQSRS QESPSEEDDG SDSQDSSRSK EDSNSTESVS SSEEEAQTKN TEVESRKLTV DAYHNKPIGD QDDNDCQDGY |
| Source | Yeast |
| Target Names | DMP1 |
| Protein Names | Recommended name: Dentin matrix acidic phosphoprotein 1 Short name= DMP-1 Short name= Dentin matrix protein 1 |
| Expression Region | 17-510 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | Dentin matrix acidic phosphoprotein is an extracellular matrix protein and a member of the small integrin binding ligand N-linked glycoprotein family. This protein, which is critical for proper mineralization of bone and dentin, is present in diverse cells of bone and tooth tissues. The protein contains a large number of acidic domains, multiple phosphorylation sites, a functional arg-gly-asp cell attachment sequence, and a DNA binding domain. In undifferentiated osteoblasts it is primarily a nuclear protein that regulates the expression of |



osteoblast-specific genes. During osteoblast maturation the protein becomes phosphorylated and is exported to the extracellular matrix, where it orchestrates mineralized matrix formation. Mutations in the gene are known to cause autosomal recessive hypophosphatemia, a disease that manifests as rickets and osteomalacia. The gene structure is conserved in mammals. Two transcript variants encoding different isoforms have been described for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.