



Recombinant Human Integrator complex subunit 12 (INTS12)

Product Code	CSB-BP011759HU
Abbreviation	INTS12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CB8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAATVNLELD PIFLKALGFL HSKSKDSA EK LKALLDESLA RGIDSSYRPS QKDVEPPKIS STKNISIKQE PKISSSLPSG NNNGKVL TTE KVKKEAEKRP ADKMKSDITE GVDIPKKPRL EKPETQSSPI TVQSSKDLPM ADLSSFEETS ADDFAMEMGL ACVVCRQMMV ASGNQLVECQ ECHNLYHRDC HKPQVTDKEA NDPRLVWYCA RCTRQMKRMA QKTQKPPQKP APAVVSVTPA VKDPLVKKPE TKLKQETTFL AFKRTEVKTS TVISGNSSSA SVSSSVTSGL TGWAAFAAKT SSAGPSTAKL SSTTQNNTGK PATSSANQKP VGLTGLATSS KGGIGSKIGS NNSTTPTVPL KPPPPLTLGK TGLSRSVSCD NVSKVGLPSP SSLVPGSSSQ LSGNGNSGTS GPSGSTTSKT TSESSSSPSA SLKGPTSQES QLNAMKRLQM VKKKAAQKKL KK
Source	Baculovirus
Target Names	INTS12
Protein Names	Recommended name: Integrator complex subunit 12 Short name= Int12 Alternative name(s): PHD finger protein 22
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.