



# Recombinant Human Integrator complex subunit 12 (INTS12)

<b>Product Code</b>	CSB-MP011759HU
<b>Abbreviation</b>	INTS12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96CB8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAATVNLELD PIFLKALGFL HSKSKDSA EK LKALLDESLA RGIDSSYRPS QKDVEPPKIS STKNISIKQE PKISSSLPSG NNNGKVL TTE KVKKEAEKRP ADKMKSDITE GVDIPKKPRL EKPETQSSPI TVQSSKDLPM ADLSSFEETS ADDFAMEMGL ACVVCRQMMV ASGNQLVECQ ECHNLYHRDC HKPQVTDKEA NDPRLVWYCA RCTRQMKRMA QKTQKPPQKP APAVVSVTPA VKDPLVKKPE TKLKQETTFL AFKRTEVKTS TVISGNSSSA SVSSSVTSGL TGWAAFAAKT SSAGPSTAKL SSTTQNNTGK PATSSANQKP VGLTGLATSS KGGIGSKIGS NNSTTPTVPL KPPPPLTLGK TGLSRVSCD NVSKVGLPSP SSLVPGSSSQ LSGNGNSGTS GPSGSTTSKT TSESSSSPSA SLKGPTSQES QLNAMKRLQM VKKKAQKKL KK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	INTS12
<b>Protein Names</b>	Recommended name: Integrator complex subunit 12 Short name= Int12 Alternative name(s): PHD finger protein 22
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.