



Recombinant Human Carboxypeptidase A6 (CPA6)

Product Code	CSB-BP005880HU
Abbreviation	CPA6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4T0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	S LSGYNYEVYH SLEEIQNWMH HLNKTHSGLI HMFSIGRSYE GRSLFILKLG RRSRLKRAVW IDCGIHAREW IGPAFCQWFV KEALLTYKSD PAMRKMLNHL YFYIMPVFNV DGYHFSWTND RFWRKTRSRN SRFRCRGVDA NRNWKVKWCD EGASMHPCDD TYCGPFPESE PEVKAVANFL RKHRKHIRAY LSFHAYAQML LYPYSYKYAT IPNFRCVESA AYKAVNALQS VYGVRYRYGP ASTTLYVSSG SSMDWAYKNG IPYAFAFELR DTGYFGFLLP EMLIKPTCTE TMLAVKNITM HLLKKCP
Source	Baculovirus
Target Names	CPA6
Protein Names	Recommended name: Carboxypeptidase A6 EC= 3.4.17.1 Alternative name(s): Carboxypeptidase B
Expression Region	130-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.