



# Recombinant *Xenopus laevis* Mitogen-activated protein kinase 8 (mapk8)

<b>Product Code</b>	CSB-BP013466XBE
<b>Abbreviation</b>	mapk8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8QHK8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRSKRDSNF SVFEIGDSTF TVLKRYQNLK PIGSGAQQIV CAAFDVLER HVAIKKLSRP FQNQTHAKRA YRELVLKCV NHKNIIGLLN VFTPQKSLEE FQDLYIVMEL MDANLCQVIQ MELDHERMSY LLYQMLCGIK HLHSAGIIHR DLKPSNIVVK SDCTLKILDF GLARTAGTSF MTPYVVTRY YRAPEVILGM GYKENVDIWS VGCILGEMIK GVLFPDGDH IDQWNKVIEQ LGTPCTEFMK KLQPTVRTYV ENRPKYAGYS FEKLFDPVLF PADSEHNKLN ASQARDLLSK MLVIDASKRI SVDDALQHPY INVWYDPLEA EAPPPKIPDK QLDEREHTIE EWKELIYKEV LDWEERAKNG VIRGQPAPLG AAVTDGSAH TSSSSGDASS MSTDPTLPSD TDSSLETSAG TLGCCR
<b>Source</b>	Baculovirus
<b>Target Names</b>	mapk8
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase 8 Short name= MAP kinase 8 Short name= MAPK 8 EC= 2.7.11.24 Alternative name(s): Stress-activated protein kinase JNK1
<b>Expression Region</b>	1-426
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.