



# Recombinant Human Zinc finger protein 563 (ZNF563)

<b>Product Code</b>	CSB-BP026844HU
<b>Abbreviation</b>	ZNF563
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TA94
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDAVAFEDVA VNFTQEEWAL LGPSQKNLYR YVMQETIRNL DCIRMIWEEQ NTEDQYKNPR RNLRCHMVER FSESKDSSQC GETFSLIRDS IVNNSICPGE DPCQSAECEE VIMGHLSLNS HIRVDSGHKP HEYQEYGEKP HTHKQRGKAF SYHHSFQSRG RPHTGKKRYE CKECGKTFSS RRNLRRHMVV QGGNRPYKCK LCGKAFFWPS LLRMHERTHT GEKPYECKQC SKAFFFYSSY RRHERMHTGE KPYECKQCSK ALPDSSSYIR HERTHTGEKP YTCKQCGKAF SVSSSLRRHE TTSAEKPYE CKQCGKTFHH LGSFQIHMKR HTGDRPHKCK ICGKGFDRPS LVRYHERIHT GEKPYECKQC GKTLSSHSSSF RRHMIMHTGG GPHKCKICGK AFVYPSVCQR HEKSHSGEKP YECKQCGKAL SHSSSFRRHM VMHTGDGPNK CKVCGKAFVY PSVCQRHEKT HWRETI
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF563
<b>Protein Names</b>	Recommended name: Zinc finger protein 563
<b>Expression Region</b>	1-476
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.