



Recombinant Human Zinc finger protein 707 (ZNF707)

Product Code	CSB-EP026975HU-B
Abbreviation	ZNF707
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96C28
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDMAQEPVTF RDVAIYFSRE EWACLEPSQR ALYRDVMLDN FSSVAALGFC SPRPDLVSRL EQWEEPWVED RERPEFQAVQ RGPRPGARKS ADPKRPCDHP AWAHKKTHVR RERAREGSSF RKGFRLDTDD GQLPRAAPER TDAKPTAFPC QVLTQRCGRR PGRRERRKQR AVELSFICGT CGKALSCHSR LLAHQTVHTG TKAFECPECG QTFRWASNLQ RHQKNHTREK PFCCEACGQA FSLKDRLAQH RKVHTEHPY SCGDGCKAFK QKSNLLRHQL VHTGERPFYC ADCGKAFRTK ENLSHHQRVH SGEKPYTCAE CGKSFRWPKG FSIHRRLLHT KRFYECGHCG KGFRHLGFFT RHQRTHRHGE V
Source	E.coli
Target Names	ZNF707
Protein Names	Recommended name: Zinc finger protein 707
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.