



Recombinant Human Zinc finger protein 333 (ZNF333)

Product Code	CSB-YP026681HU
Abbreviation	ZNF333
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96JL9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESVTFEDVA VEFIQEWALL DSARRSLCKY RMLDQCRTLA SRGTPPCKPS CVSQLGQRAE PKATERGILR ATGVAWESQL KPEELPSMQD LLEEASSRDM QMGPGLFLRM QLVPSIEERE TPLTREDRPA LQEPWWSLGC TGLKAAMQIQ RVVIPVPTLG HRNPWVARDS AVPARDPAWL QEDKVEEEAM APGLPTACSQ EPVTFADVAV VFTPEEWVFL DSTQRSLYRD VMLENYRNLA SVADQLCKPN ALSYLEERGE QWTTDRGVLS DTCAEPQCQP QEAIPSQDTF TEILSIDVKG EQPQPGEKLY KYNELEKPFN SIEPLFQYQR IHAGEASCEC QEIRNSFFQS AHLIVPEKIR SGDKSYACNK CEKSFYSSD LIRHEKHTHA EKCFCDCQECG QAFKYSSNLR RHMRTHTGEK PFECSSQCGKT FTRNFNLILH QRNHTGEKPY ECKDCGKAFN QPSSLRSHVR THTGEKPFEC SQCGKAFREH SSLKTHLRTH TREKPYECNQ CGKPFRTSTH LNVHKRIHTG EKLYECATCG QVLSRLSTLK SHMRHTHTGEK PYVCQECGRA FSEPSSLRKH ARTHSGKKPY ACQECGRAFG QSSHILIVHVR THSAGRPYQC NQCEKAFRHS SSLTVHKRTH VGRETIRNGS LPLSMHPYC GPLAN
Source	Yeast
Target Names	ZNF333
Protein Names	Recommended name: Zinc finger protein 333
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.