



# Recombinant Human Zinc finger protein 366 (ZNF366)

<b>Product Code</b>	CSB-MP026702HU
<b>Abbreviation</b>	ZNF366
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N895
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MQKEMKMIKD EDVHFDLAVK KTPSFPHCLQ PVASRGKAPQ RHPFPEALRG            PFSQFRYEPP PGDLDFGPGV FEGAGSRKRK SMPTKMPYNH PAEEVTLALH            SEENKNHGLP NPLLLFPQPP RPKYDSQMID LCNVGFQFYR SLEHFGGKPV            KQEPIKPSAV WPQPTPTPFL PTPYYPYPKV HPGLMFPPFV PSSSPFFFSR            HTFLPKQPPE PLLPRKAEPQ ESEETKQKVE RVDVNVQIDD SYYVDVGGSQ            KRWQCPTCEK SYTSKYNLVT HILGHSGIKP HACTHCGKLF KQLSHLHTHM            LTHQGTRPHK CQVCHKAFTQ TSHLKRHMMQ HSEVKPHNCR            VCGRGFAYPS ELKAHEAKHA SGRENICVEC GLDFPTLAQL KRHLTTHRGP            IQYNCSECDK TFQYPSQLQN HMMKHKDIRP YICSECGMEF VQPHHLKQHS            LTHKGVKEHK CGICGREFTL LANMKRHVLI HTNIRAYQCH LCYKSFVQKQ            TLKAHMIVHS DVKPFKCKLC GKEFNRMHNL MGHMHLHSDS KPFKCLYCPS            KFTLKGNLTR HMKVKHGVME RGLHSQGLGR GRIALAQTAG VLRSLEQEEP            FDLSQKRRAK VPVFQSDGES AQQSHCHEEE EEDNCYEVVEP YSPGLAPQSQ            QLCTPEDLST KSEHAPEVLE EACKEEKEDA SKGEWEKRSK GDLGAEGGQE            RDCAGRDECL SLRAFQSTRR GPSFSDYLYF KHRDESLKEL LERKMEKQAV            LLGI         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF366
<b>Protein Names</b>	Recommended name: Zinc finger protein 366
<b>Expression Region</b>	1-744
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.