



Recombinant Human Zinc finger protein 366 (ZNF366)

| | |
|--------------------------|---|
| Product Code | CSB-YP026702HU |
| Abbreviation | ZNF366 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8N895 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | <p>MQKEMKMIKD EDVHFDLAVK KTPSFPHCLQ PVASRGKAPQ RHPFPEALRG PFSQFRYEPP PGDLDFGPGV FEGAGSRKRK SMPTKMPYNH PAEEVTLALH SEENKNHGLP NLPLLFPQPP RPKYDSQMID LCNVGFQFYR SLEHFGGKPV KQEPIKPSAV WPQPTPTPFL PTPYYPYPKV HPGLMFPPFV PSSSPFFFSR HTFLPKQPPE PLLPRKAEPQ ESEETKQKVE RVDVNVQIDD SYYVDVGGSQ KRWQCPTCEK SYTSKYNLVT HILGHSGIKP HACTHCGKLF KQLSHLHTHM LTHQGTRPHK CQVCHKAFTQ TSHLKRHMMQ HSEVKPHNCR VCGRGFAYPS ELKAHEAKHA SGRENICVEC GLDFPTLAQL KRHLTTHRGP IQYNCSECDK TFQYPSQLQN HMMKHKDIRP YICSECGMEF VQPHHLKQHS LTHKGVKEHK CGICGREFTL LANMKRHVLI HTNIRAYQCH LCYKSFVQKQ TLKAHMIVHS DVKPFKCKLC GKEFNRMHNL MGHMHLHSDS KPFKCLYCPS KFTLKGNLTR HMKVKHGVME RGLHSQGLGR GRIALAQTAG VLRSLEQEEP FDLSQKRRAK VPVFQSDGES AQQSHCHEEE EEDNCYEVEP YSPGLAPQSQ QLCTPEDLST KSEHAPEVLE EACKEEKEDA SKGEWEKRSK GDLGAEGGQE RDCAGRDECL SLRAFQSTRR GPSFSDYLYF KHRDESLKEL LERKMEKQAV LLGI</p> |
| Source | Yeast |
| Target Names | ZNF366 |
| Protein Names | Recommended name: Zinc finger protein 366 |
| Expression Region | 1-744 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.