



Recombinant Human Zinc finger protein 555 (ZNF555)

Product Code	CSB-MP026835HU
Abbreviation	ZNF555
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEP9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSVVFEDVA VDFMLEEWAL LDSAQRDLR DVMLETFQNL ASVDDDETQFK ASGSVSQQDI YGEKIPKESK IATFTRNVSW ASVLGKIWDS LSIEDQTTNQ GRNLSRNHGL ERLCESNDQC GEALSQIPHL NLYKKIPPGV KQYEYNTYGK VFMHRRTSLK SPITVHTGHK PYQCQCEGQA YSCRSHLRMH VRTHNGERP VCKLCGKTFP RTSSLNRHVR IHTAEKTYEC KQCGKAFIDF SSLTSHLRSH TGEKPYKCKE CGKAFSYSST FRRHTITHTG EKPYKCKECA EAFSYSSTFR RHMISHTGEK PHKCKEAGEA FSYSSAFRRH MITHTGEKPY ECKQCGKTFI YLQSFRRHER IHTGEKPYEC KQCGKTFIYP QSFRRHERTH GGEKPYECNQ CGKAFSHPSS FRGHMRVHTG EKPYECKQCG KTFNWPISLR KHMRTHTREK PYECKQCGKA FLSACFREQ VRMHPEDKSY ECKLCGKAFY CHISLQKHMR RHTAEKLYKC KQCGKAFSWP ELLQQHVTRH TVEKPYECKE CGKVFVWPSS LPIHMLRHTG EKPQCKHCG KAFNCSSSLR RHVRIHTTEK QYKCNVGHPP ANEFMCSASE KSHQERDLIK VVNMVLP
Source	Mammalian cell
Target Names	ZNF555
Protein Names	Recommended name: Zinc finger protein 555
Expression Region	1-628
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.