



Recombinant *Xanthomonas campestris* pv. *campestris* Uracil-DNA glycosylase (ung)

Product Code	CSB-MP025641XAY
Abbreviation	ung
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8P4D7
Product Type	Recombinant Protein
Immunogen Species	<i>Xanthomonas campestris</i> pv. <i>campestris</i> (strain ATCC 33913 / DSM 3586 / NCPPB 528 / LMG 568 / P 25)
Purity	>85% (SDS-PAGE)
Sequence	MTEGEGRIQL EPSWKARVGE WLLQPQMDEL SAFLRQRKAA NARVFPPGPQ IFAAFDTPF EQVKVVVLGQ DPYHGEGQAH GLCFSVLPGV PVPPSLLNIY KEIQDDLIP RPDHGYLMPW ARQGVLLNA VLTVEQGRAG AHQNKGWEGF TDHVVELNR EREGLVLLW GSYAQSKGKV IDQARHRVFK APHPSPSAH RGFLGCKHFS KTNEHLQRRG LSPIDWSLPS RAALDSLALG G
Source	Mammalian cell
Target Names	ung
Protein Names	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.