



Recombinant Human Zinc finger protein 671 (ZNF671)

Product Code	CSB-YP026937HU
Abbreviation	ZNF671
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAW3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSPVSRDAS DALQGRKCLR PRSRRPLPA AVRAHGPM AE LTDSARGCVV FEDVFVYFSR EEWELLDDAQ RLLYHDVMLE NFALLASLGI AFSRSRAVMK LERGEWPVY DQVDMTSATE REAQRGLRPG CWHGVEDEEV SSEQSIFVAG VSEVRTLMAE LESHPCDICG PILKDTLHLA KYHGGKARQK PYLCGACGKQ FWFSTDFDQH QNQPNGGKLF PRKEGRDSVK SCRVHVPEKT LTCGKGRDF SATSGLLQHQ ASLSSMKPHK STKLVSGFLM GQRYHRCGEC GKAFTRKDTL ARHQRIHTGE RPYECNECGK FFSQSYDLFK HQTVHTGERP YECSECGKFF RQISGLIEHR RVHTGERLYQ CGKCGKFFSS KSNLIRHQEV HTGARPYVCS ECGKEFSRKH TLVLHQRTHT GERPYECSEC GKAFSQSSHL NVHWRIHSSD YECSRCGKAF SCISKLIQHQ KVHSGEKPYE CSKCGKAFTQ RPNLIRHWKV HTGERPYVCS ECGREFIRKQ TLVLHQRVHA GEKL
Source	Yeast
Target Names	ZNF671
Protein Names	Recommended name: Zinc finger protein 671
Expression Region	1-534
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.