



# Recombinant Mouse Death-associated protein kinase 2 (Dapk2)

<b>Product Code</b>	CSB-EP006500MO-B
<b>Abbreviation</b>	Dapk2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8VDF3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVQASMRSPN METFKQQKVE DFYDIGEELG SGQFAIVKKC REKSTGLEYA AKFIKKRQSR ASRRGVCREE IEREVSILRQ VLHPNIITLH DVYENRTDVV LILELVSGGE LDFLAQKES LSEEEATSFQ KQILDGVNYL HTKKIAHFDL KPENIMLLDK NIPPHIKLI DFGLAHEIED GVEFKNIFGT PEFVAPEIVN YEPLGLEADM WSIGVITYIL LSGASPFLGD TKQETLANIT AVSYDFDEEF FSQTSELAKD FIRKLLVKET RKRLTIQEAL RHPWITPVDT QQAMVRRESV VNLENFKKQY VRRRWKLSFS IVSLCNHLTR SLMKKVHLRT SEDLRNCESD TEENIARRKA LHPRRRSSTS
<b>Source</b>	E.coli
<b>Target Names</b>	Dapk2
<b>Protein Names</b>	Recommended name: Death-associated protein kinase 2 Short name= DAP kinase 2 EC= 2.7.11.1 Alternative name(s): DAP-kinase-related protein 1 Short name= DRP-1
<b>Expression Region</b>	1-370
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a protein that belongs to the serine/threonine protein kinase family. This protein contains a N-terminal protein kinase domain followed by a conserved calmodulin-binding domain with significant similarity to that of death-associated protein kinase 1 (DAPK1), a positive regulator of programmed cell death. Overexpression of this gene was shown to induce cell apoptosis. It uses multiple polyadenylation sites.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.