



Recombinant Human Tsukushin (TSKU)

Product Code	CSB-MP025140HU
Abbreviation	TSKU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUA8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TRPC FPGCQCEVET FGLFDSFSLT RVDCSGLGPH IMPVPIPLDT AHLDLSSNRL EMVNESVLAG PGYTTLAGLD LSHNLLTSIS PTAFSRLRYL ESLDLSHNGL TALPAESFTS SPLSDVNLSH NQLREVSUSA FTTHSQGRAL HVDLSHNLIIH RLVPHPTRAG LPAPTIQSLN LAWNRLHAVP NLRDLPLRYL SLDGNPLAVI GPGAFAGLGG LTHLSLASLQ RPELAPSGF RELPGLQVLD LSGNPKLNWA GAEVFSGLSS LQELDLSGTN LVPLPEALLL HLPALQSVSV GQDVRCRRLV REGTYPRRPG SSPKVALHCV DTRDSAARGP TIL
Source	Mammalian cell
Target Names	TSKU
Protein Names	Recommended name: Tsukushin Short name= Tsukushi Alternative name(s): E2-induced gene 4 protein Leucine-rich repeat-containing protein 54
Expression Region	17-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.