



Recombinant Human PHD finger protein 21A (PHF21A)

Product Code	CSB-EP017911HU-B
Abbreviation	PHF21A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BD5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MELQTLQEAL KVEIQVHQKL VAQMKQDPQN ADLKKQLHEL QAKITALSEK QKRVEQLRK NLIVKQECPD KFAQIPLPQS ENKLQTAQQQ PLQQLQQQQQ YHHHHAQQSA AASPNTASQ KVTVTASMIT TKTLPLVLKA ATATMPASVV GQRPTIAMVT AINSQKAVLS TDVQNTPVNL QTSSKVTGPG AEAQVIVAKN TVTLVQATPP QPIKVPQFIP PPRLTPRPNF LPQVRPKPVA QNNIPIAPAP PPMLAAPQLI QRPVMLTKFT PTTLPTSQNS IHPVRVNGQ TATIAKTFPM AQLTSIVIAT PGTRLAGPQT VQLSKPSLEK QTVKSHTETD EKQTESRTIT PPAAPKPKRE ENPQKLAFMV SLGLVTHDHL EEIQSKRQER KRRTTANPVY SGAVFEPERK KSAVTYLNST MHPGTRKRGR PPKYNAVLGF GALTPTSPQS SHPDSPENEK TETTFTFPAP VQPVSLPSPT STDGDIHEDF CSVCRKSGQL LMCDTCSRVI HLDCLDPPLK TIPKGMWICP RCQDQMLKKE EAIPWPGTLA IVHSYIAYKA AKEEEKQKLL KWSSDLKQER EQLEQVKVQL SNSISKCMEM KNTILARQKE MHSSLEKVKQ LIRLIHGIDL SKPVDSEATV GAISNGPDCT PPANAATSTP APSPSSQSCT ANCNQGEETK
Source	E.coli
Target Names	PHF21A
Protein Names	Recommended name: PHD finger protein 21A Alternative name(s): BHC80a BRAF35-HDAC complex protein BHC80
Expression Region	1-680
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.