



Recombinant Human Zinc finger protein 689 (ZNF689)

Product Code	CSB-MP026953HU
Abbreviation	ZNF689
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CS4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAPPSAPLPA QGPGKARPSR KRGRPRALK FVDVAVYFSP EEWGCLRPAQ RALYRDVMRE TYGHLGALGC AGPKPALISW LERNTDDWEP AALDPQEYPR GLTVQRKSRT RKKNGEKEVF PPKEAPRKGK RGRRPSKPRL IPRQTSGGPI CPDCGCTFPD HQALESHKCA QNLKKPYPCP DCGRRFSYPS LLVSHRRAHS GECPYVCDQC GKRFSSQRKLN SQHQVIHTGE KPYHCPDCGR CFRRSRSLAN HRTTHTGEKP HQCPSCGRRF AYPSSLAIHQ RTHTGEKPYT CLECNRRFRQ RTALVIHQRI HTGEKPYPCP DCERRFSSSS RLVSHRRVHS GERPYACEHC EARFSQRSTL LQHQLLHTGE KPYPCPDCGR AFRRSGSLAI HRSTHTEEKL HACDDCGRRF AYPSSLASHR RVHSGERPYA CDLCSKRFAQ WSHLAQHQLL HTGEKPFPCLECGRCFRQRW SLAVHKCSPK APNCSPRSAI GGSSQRGNAH
Source	Mammalian cell
Target Names	ZNF689
Protein Names	Recommended name: Zinc finger protein 689
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.