



Recombinant Human Zinc finger protein 668 (ZNF668)

Product Code	CSB-EP026933HU-B
Abbreviation	ZNF668
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96K58
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVEAAEARS PAPGYKRSGR RYKCVSCTKT FPNAPRAARH AATHGPADCS EEVAEVKPKP ETEAKAEEAS GEKVS GSAK PRPYACPLCP KAYKTAPELR SHGRSHTGEK PFPCPECGRR FMQPVCLRVH LASHAGELPF RCAHCPKAYG ALSKLKIHRQ GHTGERPYAC ADCGKSFADP SVFRKHRRTH AGLRPYSCER CGKAYAELKD LRNHRSHTG ERPF LCSECG KSF SRSSSLT CHQRIHAAQK PYRCPACGKG FTQLSSYQSH ERTHSGEKPF LCPRCGRMFS DPSSFRRHRQR AHEGVKPYHC EKC GKDFRQP ADLAMHRRVH TGDRPFKCLQ CDKTFVASWD LKRHALVHSG QRPFRCEECG RAFAERASLT KHSRVHSGER PFHCNACGKS FVSSSLRKH ERTHRSSEAA GVPPAQELVV GLALPVG VAG ESSAAPAAGA GLGDPPAGLL GLPPESGGVM ATQWQVVGMT VEHVECQDAG VREAPGPLEG AGEAGGEEAD EKPPQFVCRE CKETFSTMTL LRRHERSHPE LRPFPCTQCG KSFSDRAGLR KHSRTHSSVR PYTCPHCPKA FLSASDLRKH ERTHPVPMGT PTPLEPLVAL LGMPEEGPA
Source	E.coli
Target Names	ZNF668
Protein Names	Recommended name: Zinc finger protein 668
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.