



Recombinant Human Zinc finger protein 560 (ZNF560)

Product Code	CSB-MP026841HU
Abbreviation	ZNF560
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MR9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAYCLTNCYQ YSVTFEDTAV DFTQEEWILL DPVQRNLYRD VMLENYENVA KVGFLFKPS VISWLEEEEL RTLQQGVLQD WAIKHQTSVS ALQQEFWKIQ TSNGIQMDLV TFDSVAVEFT QEEWTLDDPA QRNLYSDVML ENYKNLSSVG YQLFKPSLIS WLEEEELST LPRVLQEWKM CLKTKGPALW QDNFCLKTLN GIQLARNQNG EELYDCKQCE DVFCKHPCLK TNMSTQNRGN TSECIQYAKD LLSLYNKTST IRKVSVFSKH GKSFRLLILNV QVQRKCTQDK SFEGTDYGKA FIYQSYLEAH RKTQSGEKLN EWKQCGEAFT HSTSHAVNVE THIIKNPYEC KECGKDFRYP THLNNHMQTH IGIKPYKCKH CGKTFTVPSG FLEHVRHTGT EKPYGCKEKG KAFGTSAGLI EHIRCHAREK TFKCDHCGKA FISYPSLFGH LRVHNGEKPY EHKEYGKAFG TSSGVIEDRR SNTGQKRFDG DQCGKVFVVF SSLFAHLRTH TGEKPFKCYK CGKPFTSSAC LRIHMRTHTE ERLYQCKKCG KAFTKCSYLT KHLRTHAGEK PYECMKCGKA FTERSYLTKH LRRHSGEKPY ECKKCGKAFT ERSDLTKHLR RHTGDKPYEY KDCGKAFFVS SSLVDHLRTH TGYKPYKCNA CEKAYSRSVCV LTQHLKTHAA EKTSECNACG NSFRNSMCFH DRLKTLTKIK PYKCKDCGKA FTCHSDLTNH VRIHTGEKPY KCKEKGKAFR TSSGRIQHLR THMGEKPFEC DQCGKAFASF SARIAHLKTH
Source	Mammalian cell
Target Names	ZNF560
Protein Names	Recommended name: Zinc finger protein 560
Expression Region	1-790
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.