



# Recombinant Human Zinc finger protein 560 (ZNF560)

<b>Product Code</b>	CSB-YP026841HU
<b>Abbreviation</b>	ZNF560
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96MR9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAYCLTNCYQ YSVTFEDTAV DFTQEEWILL DPVQRNLYRD VMLENYENVA KVGFLFKPS VISWLEEEEL RTLQQGVLQD WAIKHQTSVS ALQQEFWKIQ TSNGIQMDLV TFDSVAVEFT QEEWTLDDPA QRNLYSDVML ENYKNLSSVG YQLFKPSLIS WLEEEEEELST LPRVLQEWKM CLKTKGPALW QDNFCLKTLN GIQLARNQNG EELYDCKQCE DVFCKHPCLK TNMSTQNRGN TSECIQYAKD LLSLYNKTST IRKVSVFSKH GKSFRLLILNV QVQRKCTQDK SFEGTDYGKA FIYQSYLEAH RKTQSGEKLN EWKQCGEAFT HSTSHAVNVE THIIKNPYEC KECGKDFRYP THLNNHMQTH IGIKPYKCKH CGKTFTVPSG FLEHVRHTGT EKPYGCKEKG KAFGTSAGLI EHIRCHAREK TFKCDHCGKA FISYPSLFGH LRVHNGEKPY EHKEYGKAFG TSSGVIEDRR SNTGQKRFDQ DQCGKVFVVF SSSLFAHLRTH TGEKPFKCYK CGKPFTSSAC LRIHMRTHTE ERLYQCKKCG KAFTKCSYLT KHLRTHAGEK PYECMKCGKA FTERSYLTKH LRRHSGEKPY ECKKCGKAFT ERSDLTKHLR RHTGDKPYEY KDCGKAFFVS SSLVDHLRTH TGYKPYKCNA CEKAYSRSVCV LTQHLKTHAA EKTSECNACG NSFRNSMCFH DRLKTLTKIK PYKCKDCGKA FTCHSDLTNH VRIHTGEKPY KCKEKGKAFR TSSGRIQHLR THMGEKPFEC DQCGKAFASF SARIAHLKTH
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF560
<b>Protein Names</b>	Recommended name: Zinc finger protein 560
<b>Expression Region</b>	1-790
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.