



Recombinant Human Hyaluronan and proteoglycan link protein 3 (HAPLN3)

Product Code	CSB-MP010132HU
Abbreviation	HAPLN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LPF YNGFYYSNSA NDQNLGNGHG KDLLNGVKLV VETPEETLFT YQGASVILPC RYRYEPALVS PRRVRVKWWK LSENGAPEKD VLVAIGLRHR SFGDYQGRVH LRQDKEHDVS LEIQDLRLED YGRYRCEVID GLEDESLVE LELRGVVFPY QSPNGRYQFN FHEGQQVCAE QAAVVASFEQ LFRAWEEGLD WCNAGWLQDA TVQYPIMLPR QPCGGPGLAP GVRSYGPRHR RLHRYDVFCF ATALKGRVYY LEHPEKLTLT EAREACQEDD ATIAKVGQLF AAWKFHGLDR CDAGWLADGS VRYPVVHPPH NCGPPEPGVR SFGFPDPQSR LYGVYCYRQH
Source	Mammalian cell
Target Names	HAPLN3
Protein Names	Recommended name: Hyaluronan and proteoglycan link protein 3
Expression Region	18-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene belongs to the hyaluronan and proteoglycan binding link protein gene family. This protein may function in hyaluronic acid binding and cell adhesion.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.