



# Recombinant Mouse Phospholipase A1 member A (Pla1a)

<b>Product Code</b>	CSB-BP018084MO
<b>Abbreviation</b>	Pla1a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8VI78
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	NVPPT TQPKCTDFQS ANLLRGTNLK VQFLLFTPSD PSCGQLVEEG SDIRSSEFNA SLGTKVIIHG FRALGTKPSW IDKFISAVLR AADANVIAVD WVYGSTGVYY SAVENVVKLS LEISRFLSKL LELGVSESSI HIIGVSLGAH VGGMVGHFYK GQLGQITGLD PAGPEYTRAS LEERLDAGDA LFVEAIHTDT DNLGIRIPVG HVDYFVNGGQ DQPGCPAFFH AGYNYLICDH MRAVHLYISA LENTCPLMAF PCASYKAFLA GDCLDCFNPF LLSCPRIGLV ERGGVMIEPL PKEVKVYLLT TSSAPYCVHH SLVEFYLKEK RKKDTSIEVT FLSNNVTSSV KITIPKQQL E GRGVMAHPNP QCQINQVKLK FQVSSRVWRK DRTPVVGTF C TAPLPVNDSK KTVCIPEPVR LQAGVPAFQD LKIACV
<b>Source</b>	Baculovirus
<b>Target Names</b>	Pla1a
<b>Protein Names</b>	Recommended name: Phospholipase A1 member A EC= 3.1.1.- Alternative name(s): Phosphatidylserine-specific phospholipase A1 Short name= PS-PLA1
<b>Expression Region</b>	26-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.