



# Recombinant Human Zinc finger protein 606 (ZNF606)

<b>Product Code</b>	CSB-YP026886HU
<b>Abbreviation</b>	ZNF606
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WXB4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAINPWASW GALTDQSWGM TAVDPWASWA LCPQYPAWHV EGSLEEGRRA TGLPAAQVQE PVTFKDVAVD FTQEEWGQLD LVQRTLYRDV MLETYGHLLS VGNQIAKPEV ISLLEQGEEP WSVEQACPQR TCPEWVRNLE SKALIPAQSI FEEEQSHGMK LERYIWDDPW FSRLEVLGCK DQLEMYHMNQ STAMRQMVFM QKQVLSQRSS EFCGLGAEFS QNLNFVPSQR VSQIEHFYKP DTHAQSWRCD SAIMYADKVT CENNDYDKTV YQSIQPIYPA RIQTGDNLFK CTDAVKSFNH IIHFGDHKGI HTGEKLYEYK ECHQIFNQSP SFNEHPRLHV GENQYNYKEY ENIFYFSSFH EHQKIGTVEK AYKYNEWWEKV FGYSFSLTQH TSTYTAEKPY DYNECGTSFI WSSYLIQHKK THTGEKPYEC DKCGKVFRNR SALTKHERTH TGIKPYECNK CGKAFSWNSH LIVHKRIHTG EKPVCNECG KSFNWNSHLI GHQRTHTGEK PFECTECGKS FSWSSHIAH MRMHTGEKPF KCDECEKAFR DYSALSKHER THSGAKPYKC TECGKSFSWS SHLIAHQRTH TGEKPYNCQE CGKAFRERSA LTKHEIHSK IKPYECNKCG KSCSQMAHLV RHQRTHTGEK PYECNKCGKS FSQSCHLVAH RRIHTGEKPY KCNQCERSFN CSSHLIAHRR THTGEKPYRC NECGKAFNES SSLIVHLRNH TGEKPYKCNH CEKAFCKNSS LIIHQRMHSG EKRFICSECG KAFSGHSALL QHQRNHSEEK LN
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF606
<b>Protein Names</b>	Recommended name: Zinc finger protein 606 Alternative name(s): Zinc finger protein 328
<b>Expression Region</b>	1-792
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.