



Recombinant *Macaca fascicularis* Deleted in azoospermia-like (DAZL)

Product Code	CSB-EP006514MOV-B
Abbreviation	DAZL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95192
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	MSAANPETPN STISREASTQ SSSAAASQGY VLPEGKIMPN TVFVGGIDVR MDETEIRSFY ARYGSVKEVK IITDRTGVSK GYGFVSFFND VDVQKIVESQ INFHGKCLKL GPAIRKQNLK AYHVQPRPLV FNHPPPPQFQ NVWTNPNTET YMQPPTTMNP ITQYVQAYPT YPNSPVQVIT GYQLPVYNYQ MPPQWPVGEQ RSYVVPPAYS AVNYHCNEVD PGAEVVPNEC SVHEATPPSG NGPQKKSVDK SIQTVVSLCF NPENRLRNTV VTQDDYFKDK RVHHRRSRA MLKSV
Source	<i>E.coli</i>
Target Names	DAZL
Protein Names	Recommended name: Deleted in azoospermia-like Alternative name(s): DAZ-like autosomal Deleted in azoospermia-like 1
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.