



Recombinant Human Zinc finger protein 583 (ZNF583)

Product Code	CSB-MP026865HU
Abbreviation	ZNF583
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96ND8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSKDLVTFGD VAVNFSQEEW EWLNPAQRNL YRKVMLENYR SLVSLGVSVS KPDVISLLEQ GKEPWMVKKE GTRGPCPDWE YVFNSEFSS KQETYEESK VVTVGARHLS YSLDYPRLRE DCQSEDWYKN QLGSQEVHLS QLIITHKEIL PEVQNKEYNK SWQTFHQDTI FDIQQSFPTK EKAHKHEPQK KSYRKKSVEM KHRKVYVEKK LLKCNDCEKV FNQSSSLTLH QRIHTGEKPY ACVECGKTFS QSANLAQHKR IHTGEKPYEC KECRKAFSQN AHLAQHQVRH TGEKPYQCKE CKKAFSQIAH LTQHQRVHTG ERPFECIECG KAFSNGSFLA QHQRIHTGEK PYVCNVCCKA FSHRGYLIVH QRIHTGERPY ECKEKRKAFS QYAHLAQHQR VHTGEKPYEC KVCRKAFSQI AYLDQHQVRH TGEKPYECIE CGKAFSNSSS LAQHQRSHGT EKPYMCKEKR KTFQSNAGLA QHQRIHTGEK PYECNVCCKA FSYSGSLTLH QRIHTGERPY ECKDCRKSFR QRAHLAHHHR IHTMESFLTL SSPSPSTSNQ LPRPVGFIS
Source	Mammalian cell
Target Names	ZNF583
Protein Names	Recommended name: Zinc finger protein 583 Alternative name(s): Zinc finger protein L3-5
Expression Region	1-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.