



Recombinant Human RING finger protein 169 (RNF169)

Product Code	CSB-MP019855HU
Abbreviation	RNF169
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NCN4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAGPSTRA SSAAAAAALS RRGRRGRGRCDE TAAAKTGAPG PASGPSLLVL SPPLLQPPLP PRPEESGCAG CLEPPGEAAA LPCGHSLCRG CAQRAADAAG PGCPRCRARG PGWARRRRARD DGQADSEVLG ECARRSQPER CRPRRDGGAA AAGPRPEQEP RAAPAEPDFI FRAPIKLSKP GELREEYESL RKLREEKLQE EKPSEDQIHK LLPEDTETGK RKMDEQKKRD EPLVLKTNLE RCPARLSDSE NEEPSRGQMT QTHRSAFVSK NNSYSLAFLA GKLNSKVERS QSCSDTAQER AKSRVRAVPG NKAKVTTMTP ASNPIIGVLL STQNNRCVSA PDLTIEKRLP FSSLSSLASL HKPERSVSPE SNDSISEELN HFKPIVCSPC TPPKRLPDGR VLSPLIKST PRNLNRSQK QTSYEASPRI LKKWEQIFQE RQIKKTLKA TLTSLAPEMG EELLGSEGIH SSKEKPLVAV NTRLGGQVL SEYTGPTSAD LDHFPSVSQT KAEQSDSNKS STEIPLETCC SSELKGGGSG TSLEREQFEG LGSTPDAKLD KTCISRAMKI TTVNSVLPQN SVLGGVLKTK QQLKTLNHFD LTNGVLVESL SEEPLPSLRR GRKRHCKTKH LEQNGSLKKL RQTSGEVGLA PTDPVLREME QKLQEQEEDR QLALQLQRMF DNERRTVSRR KGSVDQYLLR SSNMAGAK
Source	Mammalian cell
Target Names	RNF169
Protein Names	Recommended name: RING finger protein 169
Expression Region	1-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.