



Recombinant Human Zinc finger protein 280C (ZNF280C)

Product Code	CSB-MP026635HU
Abbreviation	ZNF280C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8ND82
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDDDKPFQPK NISKMAELFM ECEEELEPW QKKVEETQDE DDELIFVGE ISSSKPAISN ILNRGHSSSS SKGIKSEPHS PGIPEIFRTA SQRCRDPPSN PVAASPRFHL VSKSSQSSVT VENASKPDFT KNSQVGS DNS SILLFDSTQE SLPPSQDIPA IFREGMKNTS YVLKHPSTSK VNSVTPKKPK TSEDVPQINP STSLPLIGSP PVTSSQVMLS KGTNTSSPYD AGADYLRACP KCN VQFNLLD PLKYHMKHCC PDMITKFLGV IVKSERPCDE DKTDSETGKL IMLVNEFYYG RHEGVTEKEP KTYTTFKCFS CSKVLKNNIR FMNHMKHHLE LEKQNNESWE NHTTCQH CYR QYPTPFQLQC HIESTHTPHE FSTICKICEL SFETE HILLQ HMKDTHKPGE MPYVCQVCQF RSSTFSDVEA HFRAAHENTK NLLCPFCLKV SKMATPYMNH YMKHQKKG V H RCPKCR LQFL TSKEKA E HKA QHRTFIKPK E LEGLPPGAKV TIRASLGPLQ SKLPTAPFGC APGTSFLQVT PPTSQNTTAR NPRKSNASRS KTSKLHATTS TASKVNTSKP RGRIAKSKAK PSYKQKRQRN RKNKMSLALK NIRCRRGIHK CIECHSKIKD FASHFSIYIH CSFCKYNTNC NKAFVNHMMS SHSNHPGKRF CIFKKHSGTL RGITLVCLKC DFLADSSGLD RMAKHLSQRK THTCQVIEN VSKSTSTSEP TTGCSLK
Source	Mammalian cell
Target Names	ZNF280C
Protein Names	Recommended name: Zinc finger protein 280C Alternative name(s): Suppressor of hairy wing homolog 3 Zinc finger protein 633
Expression Region	1-737
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.