



Recombinant Human Paraspeckle component 1 (PSPC1)

Product Code	CSB-MP018937HU
Abbreviation	PSPC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMLRGNLQKV RIEKNPARLR ALESVGESE PAAAAAMALA LAGEPAPPAP APPEDHPDEE MGFTIDIKSF LKPGEKTYTQ RCRLFVGNLP TDITEEDFKR LFERYGEPSE VFINRDRGFG FIRLESRTLA EIAKAELDGT ILKSRPLRIR FATHGAALTV KNLSPVVSNE LLEQAFSQFG PVEKAVVVVD DRGRATGKGF VEFAAKPPAR KALERCGDGA FLLTTTPRPV IVEPMEQFDD EDGLPEKLMQ KTQQYHKERE QPPRFAQPGT FEFYASRWK ALDEMEKQQR EQVDRNIREA KEKLEAEMEA ARHEHQLMLM RQDLMRRQEE LRRLEELRNQ ELQKRKQIQL RHEEEHRRRE EEMIRHREQE ELRRQQEGFK PNYMENREQE MRMGDMGPRG AINMGDAFSP APAGNQGPPP MMGMNMNNRA TIPGPPMGPG PAMGPEGAAN MGTPMMPDNG AVHNDRFPQG PPSQMGSPMG SRTGSETPQA PMSGVGPVSG GPGGGFGRGSQ GGNFEGPNKR RRY
Source	Mammalian cell
Target Names	PSPC1
Protein Names	Recommended name: Paraspeckle component 1 Alternative name(s): Paraspeckle protein 1
Expression Region	1-523
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.