



Recombinant Chicken Interferon regulatory factor 3 (IRF3)

Product Code	CSB-EP011818CH-B
Abbreviation	IRF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90643
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MAALDSEGDA QKLRFGPWLL NAVSSGLYRG LCWIDPDRRI FRIPWKHNAR KDV TSSDVEI FKAWAKASGR YEGNAEDPAK WKTNFRCALR STHMFMLLED RSKCNDDPHK VYAVASGVPN DRGSGGPVAG ALQQQPQLLL NHHDLALENT PTDSTEGVAA AALTQVDLDL LQSVLQHCNI SALGSQPTLW AHTGDALPED ALLLPGQDGC LPGPQFQDWR QLEEPLLLGN QPLTGGGCGQ DGAGALPVSE ECAIPAPSPA EELLFQSANP APPPPAGDIG GLPPLLDITI YYRGKMVYQE QVDDSRCVLA YQPLDPAVAE QRLVLFPSA SLPDPRQRRY TEDLLEVAGL RLEQRAGQLL ATRLKCKVF WALSQLEGG EPPLNLLHRD QETTIFDFRV FCTELRDFRD SRRERSPDFT IFLCFGQCFS STKPKEKSLI LVKLVQFCE YWYEQVQRGG ASSLN SGNVS LQLSDSFNLF ELIEQYHMQT D
Source	E.coli
Target Names	IRF3
Protein Names	Recommended name: Interferon regulatory factor 3 Short name= IRF-3
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	IRF3 encodes interferon regulatory factor 3, a member of the interferon regulatory transcription factor (IRF) family. IRF3 is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other interferon-induced genes.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.