



Recombinant Human C1GALT1-specific chaperone 1 (C1GALT1C1), partial

| | |
|--------------------------|---|
| Product Code | CSB-BP003490HU |
| Abbreviation | C1GALT1C1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q96EU7 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Source | Baculovirus |
| Target Names | C1GALT1C1 |
| Protein Names | Recommended name: C1GALT1-specific chaperone 1 Alternative name(s): C38H2-like protein 1 Short name= C38H2-L1 Core 1 beta1,3-galactosyltransferase 2 Short name= C1Gal-T2 Short name= C1GalT2 Short name= Core 1 beta3- |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |