



# Recombinant Human Zinc finger protein 300 (ZNF300)

<b>Product Code</b>	CSB-BP026653HU
<b>Abbreviation</b>	ZNF300
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96RE9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMKSQGLVSF KDVAVDFTQE EWQQLDPSQR TLYRDVMLEN YSHLVSMGYP VSKPDVISKL EQGEEPWIIK GDISNWIYPD EYQADGRQDR KSNLHNSQSC ILGTVSFHHK ILKGVTRDGS LCSILKVCQG DGQLQRFLN QDKLFRQVTF VNSKTVTEAS GHKYNPLGKI FQECIETDIS IQRFHKYDAF KKNLKPNI DL PSCYKSNSRK KPDQSFGGGK SSSQSEPNSN LEKIHNGVIP FDDNQCGNVF RNTQSLIQYQ NVETKEKSCV CVTCGKAKAFK KSQLIVHQRI HTGKKPYDCG ACGKAFSEKF HLVVHQRTHHT GEKPYDCSEC GKAFSQKSSL IIHQRVHTGE KPYECSECGK AFSQKSPLII HQRIHTGEKP YECRECGKAF SQKSQLIHH RAHTGEKPYE CTECGKAFCE KSHLIIHKRI HTGEKPYKCA QCEEAFSRKT ELITHQLVHT GEKPYECTEC GKTF SRKSQL IIHQRTHTGE KPYKCSECGK AFCQKSHLIG HQRIHTGEKP YICTECGKAF SQKSHLPGHQ RIHTGEKPYI CAECGKA FSQ KSDLVLHQRI HTGERPYQCA ICGKAFIQKS QLTVHQRIHT VVKS
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF300
<b>Protein Names</b>	Recommended name: Zinc finger protein 300
<b>Expression Region</b>	1-604
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.