



Recombinant Human Zinc finger protein 611 (ZNF611)

Product Code	CSB-BP026891HU
Abbreviation	ZNF611
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N823
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLREEAAQKR KGKEPGMALP QGRLTFRDVA IEFSLAEWKC LNPSQRALYR EVMLENYRNL EAVDISSKCM MKEVLSTGQG NTEVIHTGTL QRHESHHIGD FCFQEIEKEI HDIEFQCQED ERNGLEAPMT KIKKLTGSTD QHDHRHAGNK PIKDQLGSSF YSHLPELHIF QIKGEIGNQL EKSTNDAPSV STFQRISCRP QTQISNNYGN NPLNSSLLPQ KQEVHMKREKS FQCNKSGKAF NCSSLLRKHQ IPHLGDKQYK CDVCGKLFNH EQYLACHDRC HTVEKPYKCK ECGKTFQSQES SLTCHRRRLHT GVKRYNCNEC GKIFGQNSAL LIDKAIDTGE NPYKCNECDK AFNQQSQLSH HRIHTGEKPY KCEECDKVFS RKSTIETHKR IHTGEKPYRC KVCDTAFTHW SQLARHRRHIH TAKKTYKCNE CGKTFSHKSS LVCHHRLHGG EKSYKCKVCD KAFVWSSQLA KHTRIDCGEK PYKCNECGKT FGQNSDLLIH KSIHTGEQPY KCDECEKVFS RKSSLETHKI GHTGEKPYKC KVCDKAFACH SYLAKHTRIH SGEKPYKCNE CSKTFSHRSY LVCHHRVHSG EKPYKCNECS KTFSSRSSLH CHRRLHSGEK PYKCNECGNT FRHCSSLIYH RRLHTGEKSY KCTICDKAFV RNSLLSRHTR IHTAEKPYKC NECGKAFNQQ SHLSRHHRIH TGEKP
Source	Baculovirus
Target Names	ZNF611
Protein Names	Recommended name: Zinc finger protein 611
Expression Region	1-705
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.