



Recombinant Human Zinc finger protein 565 (ZNF565)

Product Code	CSB-MP026846HU
Abbreviation	ZNF565
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9K5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAQGLVTFRD VAIEFSLEEW KCLEPAQRDL YREVTLENFG HLASLGLSIS KPDVVSLLAQ GKEPWWMIAND VTGPWCPLD SRCEKFLQKD IFEIGAFNWE IMESLKCSDL EGSDFRADWE CEGQFERQVN EECYFKQVNV TYGHMPVFQH HTSHTVRQSR ETGEKLMECH ECGKAFSRGS HLIQHQQIHT GEKPFQCKEC GKAFSRASHL VQHQRHTGE KPYDCKDCGK AFGRTSELIL HQLRHTGVKP YECKEKGKTF RQHSQILHQ RTHTEKPYV CKDCGKAFIR GSQLTVHRR HTGARPYECK ECGKAFRQHS QLTVHQRIHT GEKPYECKEC GKGFIHSSEV TRHQRIHSGE KPYECKEKGK AFRQHAQLTR HQRVHTGDRP YECKDCGKAF SRSSYLQHQ RIHTGDKPYE CKECGKAFIR VSQ LTHHQRI HTCEKPYECR ECGMAFIRSS QLTEHQRIHP GIKPYECREC GQAFILGSQL IEHYRIHTG
Source	Mammalian cell
Target Names	ZNF565
Protein Names	Recommended name: Zinc finger protein 565
Expression Region	1-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.