



# Recombinant Human Activating signal cointegrator 1 complex subunit 1 (ASCC1)

<b>Product Code</b>	CSB-BP002196HU
<b>Abbreviation</b>	ASCC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N9N2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEVLRPQLIR IDGRNYRKNP VQEQTQHEE DEEDFYQGSM ECADEPCDAY EVEQTPQGFR STLRAPSLLY NLIHLNTSND CGFQKITLDC QNIYTWKSRH IVGKRGDTRK KIEMETKTSI SIPKPGQDGE IVITGQHRNG VISARTRIDV LLDTFRRKQP FTHFLAFFLN EVEVQEGFLR FQEEVLAKCS MDHGVDSSIF QNPKKLHLTI GMLVLLSEEE IQQTCEMLQQ CKEEFINDIS GGKPLEVEMA GIEYMNDDPG MVDVLYAKVH MKDGSNRLQE LVDRVLERFQ ASGLIVKEWN SVKLHATVMN TLFKDPNAE GRYNLYTAEG KYIFKERESF DGRNILKSFA LLPRLEYNDA ISAHCNLCLP GSSDSPASAS QVAGITGVSD AYSQSLPGKS
<b>Source</b>	Baculovirus
<b>Target Names</b>	ASCC1
<b>Protein Names</b>	Recommended name: Activating signal cointegrator 1 complex subunit 1 Alternative name(s): ASC-1 complex subunit p50 Trip4 complex subunit p50
<b>Expression Region</b>	1-400
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.