



Recombinant Human Zinc finger protein 596 (ZNF596)

Product Code	CSB-MP026876HU
Abbreviation	ZNF596
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TC21
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPSPDSMTFE DIIVDFTQEE WALLDTSQRK LFQDVMLENI SHLVSIGKQL CKSVVLSQLE QVEKLSTQRI SLLQGREGVI KHQEIPFIQH IYQKGTSTIS TMRSHQEDP FLCNDLGEDF TQHIALTQNV ITYMRTKHFV SKKFGKIFSD WLSFNQHKKEI HTKCKSYGSH LFDYAFIQNS ALRPHSVTHT REITLECRVC GKTFSKNSNL RRHEMIHTGE KPHGCHLCGK AFTHCSDLRK HERTHTGEKP YGCHLCGKAF SKSSNLRRHE MIHTREKAQI CHLCGKAFTH CSDLRKHRT HLGDKPYGCL LCGKAFSKCS YLRQHERTHN GEKPYECHLC GKAFSHCSHL RQHERSHNGE KPHGCHLCGK AFTESSVLKR HERIHTGEKP YECHVCGKAF TESSDLRRHE RTHTGEKPYE CHLCGKAFNH SSVLRRHERT HTGEKPYECN ICGKAFNRSY NFRLHRRVHT GEKPYVCPLC GKAFSKFFNL RQHERHTTKK AMNM
Source	Mammalian cell
Target Names	ZNF596
Protein Names	Recommended name: Zinc finger protein 596
Expression Region	1-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.