



# Recombinant Human Zinc finger protein 296 (ZNF296)

<b>Product Code</b>	CSB-BP026648HU
<b>Abbreviation</b>	ZNF296
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WUU4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRRKAGSAP RRVEPAPAAAN PDDEMEMQDL VIELKPEPDA QPQQAPRLGP FSPKEVSSAG RFGGEPHHSP GPMPAGAALL ALGPRNPWTL WTPLTPNYPD RQPWTDKHPD LLTCGRCLQT FPLEAITAFM DHKKLGCQLF RGPSRGQGSE REELKALSCL RCGKQFTVAW KLLRHAQWDH GLSIYQTESE APEAPLLGLA EVAAAVSAVV GPAAEAKSPR ASGSGLTRRS PTCPVCKKTL SSFSNLKVHM RSHTGERPYA CDQCPYACAQ SSKLNRHKKT HRQVPPQSPL MADTSQEQAS AAPPEPAVHA AAPTSTLPCS GGEGAGAAAT AGVQEPGAPG SGAQAGPGGD TWGAIITEQR TDPANSQKAS PCKMPKSGGK SRGPGGSCEF CGKHFTNSSN LTVHRRSHTG ERPYTCEFCN YACAQSSKLN RHRRMHGMP GSTRFECPHC HVPFGLRATL DKHLRQKHPE AAGEA
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF296
<b>Protein Names</b>	Recommended name: Zinc finger protein 296 Short name= ZFP296 Alternative name(s): Zinc finger protein 342
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.