



# Recombinant Human Zinc finger protein 480 (ZNF480)

<b>Product Code</b>	CSB-BP026779HU
<b>Abbreviation</b>	ZNF480
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WV37
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLCDEKAQKR RKRKAKESGM ALPQGHITFR DVAIEFSQAE WKCLDPAQRA LYKDVMLENY RNLVSLGISL PDLNINSMLE QRREPWSGES EVKIAKNSDG RECIKGVNTG SSSYALGSNAE DKPIKKQLGV SFHLHLSELE LFPDERVING CNQVENFINH SSSVSCLQEM SSSVKTIPIFN RNDFDDSSFL PQEQKVHLRE KPYECNEHSK VFRVSSSLTK HQVIHTVEKP YKCNSCGKVF SRNSHLAEHC RIHTGEKPYK CNVCGKVFSY NSNFARHQRI HTREKPYECN ECGKVFSSNS YLARHQRIHA EEKPYKCNEC GKGFSHKSSL ANHWRIYTGE KPYKCDECGK AFYRIALLVR HQKIHTGEKP YKCNECGKVF IQNSHLAQHW RIHTGEKPYK CNECGKVFNQ LSNLARHRRRI HTGEKPYKCN ECGKAFSEYS GLSAHLVIHT GEKPYKCSEC GKAFRHKLSL TNHQRIHTGE RPYKCNECGK VFNRI AHLAR HRKIHTGEKP YKCNECGKAF SRISYLAQHW TIHMG
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF480
<b>Protein Names</b>	Recommended name: Zinc finger protein 480
<b>Expression Region</b>	1-535
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.