



Recombinant Human Zinc finger protein 480 (ZNF480)

Product Code	CSB-EP026779HU-B
Abbreviation	ZNF480
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WV37
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLCDEKAQKR RKRKAKESGM ALPQGHLTFR DVAIEFSQAE WKCLDPAQRA LYKDVMLENY RNLVSLGISL PDLNINSMLE QRREPWSGES EVKIAKNSDG RECIKGVNTG SSYALGSNAE DKPIKKQLGV SFHLHLSELE LFPDERVING CNQVENFINH SSSVSCLQEM SSSVKTPIFN RNDFDDSSFL PQEQKVHLRE KPYECNEHSK VFRVSSSLTK HQVIHTVEKP YKCNSCGKVF SRNSHLAEHC RIHTGEKPYK CNVCGKVFSY NSNFARHQRI HTREKPYECN ECGKVFSSNS YLARHQRIHA EEKPYKCNEC GKGFSHKSSL ANHWRIYTGE KPYKCDECGK AFYRIALLVR HQKIHTGEKP YKCNECGKVF IQNSHLAQHW RIHTGEKPYK CNECGKVFNQ LSNLARHRRRI HTGEKPYKCN ECGKAFSEYS GLSAHLVIHT GEKPYKCSEC GKAFRHKLSL TNHQRIHTGE RPYKCNECGK VFNRI AHLAR HRKIHTGEKP YKCNECGKAF SRISYLAQHW TIHMG
Source	E.coli
Target Names	ZNF480
Protein Names	Recommended name: Zinc finger protein 480
Expression Region	1-535
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.