



Recombinant Human Inositol hexakisphosphate kinase 3 (IP6K3)

Product Code	CSB-BP011772HU
Abbreviation	IP6K3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96PC2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVVQNSADAG DMRAGVQLEP FLHQVGGHMS VMKYDEHTVC KPLVSREQRF YESLPLAMKR FTPQYKGTVT VHLWKDSTGH LSLVANPVKE SQEPFKVSTE SAAVAIWQTL QQTTGNSGSD CTLAQWPHAQ LARSPKESPA KALLRSEPHL NTPAFSLVED TNGNQVERKS FNPWGLQCHQ AHLTRLCSEY PENKRHRFLL LENVVSQYTH PCVLDLKMGT RQHGD DASEE KKHMRKCA QSTSACLGVR ICGMQVYQTD KKYFLCKDKY YGRKLSVEGF RQALYQFLHN GSHLRRELLE PILHQLRALL SVIRSQSSYR FYSSSLLVIY DGQEPPERAP GSPHPHEAPQ AAHGSSPGGL TKVDIRMIDF AHTTYKGYWN EHTTYDGPD GYIFGLENLI RILQDIQEGE
Source	Baculovirus
Target Names	IP6K3
Protein Names	Recommended name: Inositol hexakisphosphate kinase 3 Short name= InsP6 kinase 3 EC= 2.7.4.21 Alternative name(s): Inositol hexaphosphate kinase 3
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.