



Recombinant Human Zinc finger protein 607 (ZNF607)

Product Code	CSB-MP026887HU
Abbreviation	ZNF607
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96SK3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSYGSITFGD VAIDFSHQEW EYLSLVQKTL YQEVMMENYD NLVSLAGHSV SKPDLITLLE QGKEPWMIVR EETRGECTDL DSRCEIISDG KMQLYRKHSC VTLHQRIHNG QKPYECKQCQ KSFSHLTELM VHQTIHTSEE PDQCEKFRKA FSHLTDLRKH QKINAREKPY ECEECGKVFS YPANLAQH GK VHVEKPYECK ECGEAFRTSR QLTVHHRFHY GEKPYECKEC GKAFSVYGR LSRHQSIHTGE KPFECKNCGK SFRLKAGLKV HQSIHTGEKP HECKECKGKAF RQFSHLVGHK RIHTGEKPYE CKECGKGFTC RYQLTMHQRI YSGEKHYECK ENGEAFSSGH QLTAPHTFES VEKPYKCEEC GKAFSVHGRL TRHQGIHSGK KPYECNCKGK SFRLNSSLKI HQNIHTGEKP YKCKECKGKAF SQRAHLAHHN RIHTGYKPF ECKECKGSFRC ASYLVIHERI HTGEKPYVCQ ECGKGFSYSH KLTIHRRVHT GEKPYECKEC GKAFSVSGQL TQHLHSIHS GK KPFECKNCKGK SFRFISVLKA HQNIHSAEKP YECKECKGKAF RHATSLIYHD RTHAGEKSYE CKECGETFSH ASHLIIHERI HTSDKPYECK RCGKAFHCAS YLVRHESVHA DGNPYMCCEEC GKAFNSSHEL SIHHRVHTGE KPFCNCKCRR SFRLRSILEV HQRIHI
Source	Mammalian cell
Target Names	ZNF607
Protein Names	Recommended name: Zinc finger protein 607
Expression Region	1-696
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.