



# Recombinant Human Zinc finger protein 610 (ZNF610)

<b>Product Code</b>	CSB-MP026890HU
<b>Abbreviation</b>	ZNF610
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N9Z0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLCDEEAQKR KAKESGMALP QGRLTFMDVA IEF SQEEWKS LDPGQRALYR DVMLENYRNL VFLGICLPDL SIISMLKQRR EPLILQSQVK IVKNTDGREC VRSVNTGRSC VLG SNAENKP IKNQLGLTLE AHLSELQLFQ AGRKIYRSNQ VEKFTNHRSS VSPLQKISSS FTTHIFNKYR NDLIDFPLLP QEEKAYIRGK SYEYECSEDG EVFRVRASLT NHQVIHTAEK PYKCTECGKV FSRNSHLVEH WRIHTGQKPY KCSECDKVFN RNSNLARHQR IHTGEKPHKC NECGKAFREC SGLTTHLVIH TGEKPYKCNE CGKNFRHKFS LTNHQRSHTA EKPYKCNECG KVFSLLSYLA RHQIIHSTEK PYKCNECGRA FHKRPGLMAH LLIHTGEKPY KCNECDKVFG RKLYLTNHQR IHTGERPYKC NACGKVFNQ NPHLSRHRKIH AGENSLRTLQ ME
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF610
<b>Protein Names</b>	Recommended name: Zinc finger protein 610
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.